


## QUALITY ASSURANCE FOR RADIOGRAPHY/FLUOROSCOPY

### Attenuation body PAVO

#### Key features:

- PMMA-plate, dimensions 300 x 300 x 30 mm, coated with an 1 mm Cu-layer
- Complete attenuation body equivalent to an attenuation of 25 mm Al
- Especially developed as attenuation body for test device FLURO for QA within digital and analogue Radiography and Fluoroscopy
- Easy-to-use: Simply put PAVO directly on top of / below FLURO (depending on the positioning of the tube) and make an X-ray image



#### PAVO attenuation body, facts / delivery scope:

- Attenuation body PAVO
- Weight 5.5 kg – dimensions 300 x 300 x 31 mm
- Order number: 10007

#### Perfect in combination with:

- FLURO Phantom – Order number: 10006

#### and, recommended for evaluation of the perpendicular beam:

- Beam alignment test tool BENO – Order number: 10011

#### Relevant standard:

German DIN 6868-4, 2007